



## Choose the right stethoscope to help you perform at your best.

3M™ Littmann® Stethoscopes feature adjustable headsets with soft-sealing eartips. Select Littmann mechanical stethoscopes are designed with tunable technology to give you the ability to hear different frequency sounds by simply adjusting the pressure applied to the chestpiece. Get the acoustic performance, high quality and comfort you need with a Littmann stethoscope.

	ACOUSTIC PERFORMANCE LEVELS								
	6	7-8			9-10		10+		
	<b>Lightweight II S.E. Stethoscope</b>	<b>Classic II Pediatric and Infant Stethoscopes</b>	<b>Master Classic II™ Stethoscope</b>	<b>Classic III™ Monitoring Stethoscope</b>	<b>Master Classic II™ Veterinary Stethoscope</b>	<b>Cardiology IV™ Diagnostic Stethoscope</b>	<b>Master Cardiology™ Stethoscope</b>	<b>Electronic Stethoscope Model 3200</b>	
	Ideal for taking blood pressure readings and making limited patient assessments.	Hard-working, dependable, clinical tools for physical assessment and monitoring of your patients, one long shift after another.			Provides outstanding acoustic sensitivity to help you perform detailed diagnostic auscultation. Useful in challenging and critical care environments.		High-performance auscultation. Hear more of what you need to hear with amplification of 24x, ambient noise reduction and optional connectivity.		
Approximate Length	28"/71cm	28"/71cm	27"/68cm	27"/69cm	32"/81 cm	22"/56cm 27"/69cm	22"/56cm 27"/68cm	27"/68cm	
Approximate Weight	118g	105g, 95g	160g	150g	173g	167g, 177g	175g, 185g	185g	
Warranty	2 year	3 year	3 year	5 year	3 year	7 year	7 year	2 year	
Diaphragm	Tunable	Non-Tunable (Floating)	Tunable	Tunable (2)	Tunable	Tunable (2)	Tunable	Electronic	
Approximate Chestpiece Diameter	Double-sided 2.0" x 1.8"/5.2 x 4.5cm 1.4" x 1.3"/3.5 x 3.2cm	Double-sided Pediatric 1.4"/3.4cm 1.1"/2.7cm Infant 1.1"/2.8cm 0.7"/1.7cm	Single-sided 1.7"/4.3cm	Double-sided Adult side 1.7"/4.3cm Pediatric side 1.4"/3.4cm	Single-sided 1.75"/4.4cm	Double-sided Adult side 1.7"/4.3cm Pediatric side 1.4"/3.4cm	Single-sided 1.9"/4.8cm	Single-sided 2"/5.1cm	
Features	Teardrop-shaped chestpiece Fits under blood pressure cuff	Specially designed for the smallest patients	Designed for easier grip and maneuvering	Convertible open bell	Designed for easier grip and maneuvering	Convertible open bell Dual-lumen tubing	Dual-lumen tubing Best acoustic response in Littmann mechanical line	Eliminates 85% (on average) of ambient noise Bluetooth® wireless technology Recording and playback	



The power to learn, in the palm of your hand.



For more information, please contact Asia Vet Medical Limited.

(852) 2335 9200

(852) 2335 0008

Sales@avmhc.com



# 3M™ Littmann® Stethoscope Ordering Information and Tube Color Guide

**Components**  
(Chestpiece/Eartube/Stem)

Lightweight II SE Stethoscope	Classic II Infant Stethoscope	Classic II Pediatric Stethoscope	Master Classic II™ Stethoscope	Classic III™ Monitoring Stethoscope	Master Classic II™ Veterinary Stethoscope	Cardiology IV™ Diagnostic Stethoscope	Master Cardiology™ Stethoscope	Electronic Model 3200 Stethoscope
-------------------------------	-------------------------------	----------------------------------	--------------------------------	-------------------------------------	---	---------------------------------------	--------------------------------	-----------------------------------

**Tube Color**

All Standard Finish	2450 Black	2114 Black	2113 Black	2144L Black	5620 Black	5627 Burgundy	1392 Navy	6151 (22") Black	6156 Plum			3200BK27 Black	
	2451 Burgundy	2114R Red	2113R Red	2146 Burgundy	5621 Gray	5633 Pearl Pink		6152 Black	6158 Raspberry	2160 Black		3200BU Burgundy	
	2452 Caribbean			2147 Navy	5622 Navy	5831 Plum		6154 Navy	6159 Rose Pink	2163 Burgundy		3200NB Navy	
		2124 Caribbean	2119 Caribbean		5623 Caribbean	5832 Lavender		6155 Hunter Green			2164 Navy		
	2454 Ceil Blue		2122 Raspberry		5630 Ceil Blue	5835 Turquoise							
						5839 Lemon-Lime					2167 Plum		
	2456 Pearl Pink									2178 Caribbean			
All Black				2141 Black	5803 Black	5868 Burgundy		6163 Black		2161 Black			
								6168 Navy					
All Brass Finish										2175 Black			
All Copper Finish					5809 Chocolate								
All Smoke					5811 Black	5871 Raspberry		6162 Black		2176 Black			
					5812 Dark Olive			6166 Plum	6171 Turquoise	2182 Dark Olive			
Rainbow/Brass/Brass	2157 Raspberry	2153 Caribbean			5806 Raspberry	5807 Caribbean							
Brass/Smoke/Smoke								6164 Black					
Rainbow/Black/Black					5870 Black			6165 Black					
Mirror/Smoke/Smoke					5862 Raspberry								
					5863 Navy								
Mirror/Standard/Mirror								6170 Burgundy					
								6177 Black					
Champagne/Smoke/Smoke					5861 Black	5864 Burgundy		6179 Black					
Champagne/Black/Champagne								6176 Burgundy	6190 Caribbean				
<b>COLORED STEM SERIES: (square inset represents stem color)</b>													
Black/Black/Color								6200 Black	6203 Black				
								6201 Black					
High-Polish Smoke/Black/Color								6202 Navy	6204 Black				
Rainbow/Black/Color								6205 Plum					
Champagne/Champagne/Color								6206 Hunter					
Smoke/Smoke/Color					5872 Turquoise	5875 Lime							
					5873 Gray								
Mirror/Standard/Color					5874 Caribbean								

**Components:**  
Chestpiece/Eartube/Stem

*Standard Finish*

*Black*

*Brass Finish*

*Copper Finish*

*Smoke*

*Rainbow*

*Mirror*

*Champagne*

*High-Polish Smoke*

Colors shown are approximate.

3M, Littmann, the L Littmann logo, Classic III, Cardiology IV, Master Classic II and Master Cardiology are trademarks of 3M. Used under license in Canada. Bluetooth is a registered trademark of Bluetooth SIG. © 3M 2019. All rights reserved. 70-2011-5877-4 Sept 2019