

HEINE® Diagnostic Set

Absolute precision is the basic requirement to ensure the creation of a high-quality product. At HEINE® this principle applies to the entire production process.

Complete set including:

NT 300 Double Charger (302550)

- Automatic charging process and overcharge protection
- Includes two NT 3.5 V Battery Handles including batteries (2 x 302535)
- With short-circuit protection and a continuous ceramic brightness control

BETA® 200 Direct Ophthalmoscope (301135)

- A spherical optics eliminate reflexes and guarantee a full view is even with small pupils

Slit Illumination Otoscope Head (300035)

- Homogeneous light and a swivel lens with 1.6x magnification

Otoscope Specula

- High quality chrome plated metal construction, which is easy to clean, disinfect and sterilise
- Metal Speculum, 4 x 57 mm (300357)
- Metal Speculum, 6 x 65 mm (300365)
- Metal Speculum, 7 x 90 mm (300390)

