

# The accurate Catalyst SDMA Test you know and trust. Now upgraded.

**Easier to use. Faster results. And simple freezer storage with other routine slides.**

Kidneys are earlier indicators of many aspects of kidney health. By including the new Catalyst® SDMA Test with every chemistry panel, it's now easier than ever to hear everything kidneys have to say.



## What's changed

- + **Convenient storage:** Now in the freezer, making it easier to grab and go when getting other Catalyst® slides.
- + **Simpler workflow and faster results:** No more reagent cups. Now that the reagent has been incorporated into the slide, you can simply add the slide to your chemistry run for faster results.

## What's staying the same

- + The industry-leading, reliable accuracy you've come to expect so you can gain deeper insights into kidney function.
- + Flexibility to add to any panel in any order.
- + Same quantity, less bulk: You'll still get 12 tests per box, but smaller packaging means easier storage and less waste.

